




Phospho-TYK2 (Tyr1054) Antibody

Product Code	CSB-PA154201
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29597
Immunogen	Peptide sequence around phosphorylation site of tyrosine 1054 (H-E-Y(p)-Y-R) derived from Human TYK2.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of TYK2 only when phosphorylated at tyrosine 1054.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography usi
Clonality	Polyclonal
Alias	tyrosine kinase 2; JTK1;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TYK2

Image



Western blot analysis of extracts from HT29 cells untreated(lane 1) or treated with Anisomycin(lane 2) using TYK2(Phospho-Tyr1054) Antibody.

Product Modify	Phospho-Tyr1054
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.