




Phospho-HDAC4 (Ser632) Antibody

Product Code	CSB-PA247261
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P56524
Immunogen	Peptide sequence around phosphorylation site of serine 632 (A-Q-S(p)-S-P) derived from Human HDAC4.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of HDAC4 only when phosphorylated at serine 632.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using
Clonality	Polyclonal
Alias	HD4
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	HDAC4

Image



Western blot analysis of extracts from 293 cells untreated(lane 1) or treated with EGF(lane 2) using HDAC4(Phospho-Ser632) Antibody.

Product Modify	Phospho-Ser632
-----------------------	----------------