





Phospho-HSP90AB1 (Ser254) Antibody

Product Code	CSB-PA060032
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08238
Immunogen	Peptide sequence around phosphorylation site of serine 254 (V-G-S(p)-D-E) derived from Human HSP90B.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of HSP90B only when phosphorylated at serine 254.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography usi
Clonality	Polyclonal
Alias	HS90B; HSP 84; HSP90-beta; HSP90AB1; HSPC2
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	HSP90AB1

Image



Western blot analysis of extracts from HeLa cells untreated (lane 1) or treated with hot shock (lane 2) using HSP90B(Phospho-Ser254) Antibody.



Western blot analysis of extracts from HeLa cells, treated with calf intestinal phosphatase (CIP), using HSP90B (Phospho-Ser254) Antibody.

Product Modify Phospho-Ser254



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.