




Phospho-PXN (Tyr31) Antibody

Product Code	CSB-PA047961
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49023
Immunogen	Peptide sequence around phosphorylation site of tyrosine 31 (T-P-Y(p)-S-Y) derived from Human Paxillin.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous level of Paxillin only when phosphorylated at tyrosine 31.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography usi
Clonality	Polyclonal
Alias	PAXI; PXN;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PXN

Image



Western blot analysis of extracts from HT29 cells untreated or treated with PMA using Paxillin(Phospho-Tyr31) Antibody.

Product Modify	Phospho-Tyr31
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.