




Phospho-KIT (Tyr721) Antibody

Product Code	CSB-PA449794
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10721
Immunogen	Peptide sequence around phosphorylation site of tyrosine 721 (N-E-Y(p)-M-D) derived from Human c-Kit.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of c-Kit only when phosphorylated at tyrosine 721.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography usi
Clonality	Polyclonal
Alias	KIT , CD117, SCFR, PBT
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	KIT

Image



Western blot analysis of extracts from HeLa cells untreated or treated with IFN- α using c-Kit(Phospho-Tyr721) Antibody.

Product Modify	Phospho-Tyr721
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.