




Phospho-GAP43 (Ser41) Antibody

Product Code	CSB-PA011221
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17677
Immunogen	Peptide sequence around phosphorylation site of Ser41 (Q-A-S(p)-F-R) derived from Human GAP43.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of GAP43 only when phosphorylated at Ser41.
Tested Applications	ELISA,IF;IF:1:100-1:200
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography
Clonality	Polyclonal
Alias	B-50; BASP2; NEUM; PP46; axonal membrane protein GAP-43
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GAP43

Image



Immunofluorescence staining of methanol-fixed HeLa cells using GAP43(Phospho-Ser41) Antibody.

Product Modify	Phospho-Ser41
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.