




Phospho-ITGB3 (Tyr785) Antibody

Product Code	CSB-PA065810
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05106
Immunogen	Peptide sequence around phosphorylation site of tyrosine 785 (I-T-Y(p)-R-G) derived from Human Integrin b3.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous level of Integrin b3 only when phosphorylated at tyrosine 785.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using
Clonality	Polyclonal
Alias	CD61 antigen; GP3A; GPIIIa; ITB3; Platelet membrane glycoprotein IIIa
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ITGB3

Image



Western blot analysis of extracts from HUVEC cells using Integrin b3(Phospho-Tyr785) Antibody and the same antibody preincubated with blocking peptide .

Product Modify	Phospho-Tyr785
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.