




Phospho-GRIA2 (Ser880) Antibody

Product Code	CSB-PA824200
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P42262
Immunogen	Peptide sequence around phosphorylation site of serine 880 (I-E-S(p)-V-K) derived from Human Glutamate receptor 2.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of Glutamate receptor 2 (Precursor) only when phosphorylated at serine 880.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using
Clonality	Polyclonal
Alias	GRIA2; GluR-2; GluR-B; GluR-K2; AMPA 2
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GRIA2

Image



Western blot analysis of extracts from mouse brain tissue using Glutamate receptor 2(Precursor)(phospho-Ser880) antibody(Lane 2) and the same antibody preincubated with blocking peptide(Lane1).

Product Modify	Phospho-Ser880
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.