



Phospho-SUFU (Ser342) Antibody

Product Code	CSB-PA229368
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UMX1
Immunogen	Peptide sequence around phosphorylation site of Serine 342 (A-P-S(p)-R-K) derived from Human Sufu.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous level of Sufu protein only when phosphorylated at serine 342.
Tested Applications	ELISA, WB; WB: 1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	SUFUH
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SUFU
Product Modify	Phospho-Ser342
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.