





NFKB1 (Ab-907) Antibody

Product Code	CSB-PA696291
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19838
Immunogen	Peptide sequence around aa.905~909 (P-L-S-P-A) derived from Human NFkB-p105.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous level of total NFkB-p105 protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:1000,IHC:1:50-1:200
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	DNA-binding factor KBF1; EBP-1; NF-kappa-B1 p84/NF-kappa-B1 p98; NFKB1; NFkB-p50
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NFKB1

Image



Western blot analysis of extracts from 293 and MCF cells using NFKB-p105(Ab-907) Antibody.



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using NFkB-p105(Ab-907) Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.