

Image








MYC (Ab-373) Antibody

Product Code	CSB-PA552586
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01106
Immunogen	Peptide sequence around aa.371~375 (K-R-S-F-F) derived from Human Myc.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total Myc protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitopespecific peptide.
Clonality	Polyclonal
Alias	c-myc
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MYC
_	



Western blot analysis of extract from JK, HL-60, HUVEC and Hela cells using Myc(Ab-373) Antibody