





STAT3 (Ab-727) Antibody

Product Code	CSB-PA088158
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P40763
Immunogen	Peptide sequence around aa. 725~729 (P-M-S-P-R) derived from Human STAT3.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total STAT3 protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:1000,IHC:1:50-1:200
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	Acute-phase response factor
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	STAT3

Image



Western blot analysis of extracts from HeLa cells using STAT3(Ab-727) Antibody.



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using STAT3(Ab-727) Antibody(left) or the same antibody preincubated with blocking peptide(right).