





STAT5A (Ab-694) Antibody

Product Code	CSB-PA094290
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P42229
Immunogen	Peptide sequence around aa. 692~696 (D-G-Y-V-K) derived from Human STAT5A.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total STAT5A protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:1000,IHC:1:50-1:200
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	MGF, MPF, Mammary gland factor, STA5A, STAT5
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	STAT5A

Image



Western blot analysis of extracts from HeLa and K562 cells using STAT5A(Ab-694) Antibody.



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using STAT5A(Ab-694) Antibody(left) or the same antibody preincubated with blocking peptide(right).

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.