





Bad (Ab-112) Antibody

Product Code	CSB-PA103572
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q61337
Immunogen	Peptide sequence around aa.110~114 (H-S-S-Y-P) derived from Mouse BAD.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total BAD protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:1000,IHC:1:50-1:200
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	Bbc2
Product Type	Polyclonal Antibody
Immunogen Species	Mus musculus (Mouse)
Target Names	Bad

Image



Western blot analysis of extracts from 293 cells using BAD(Ab-112) Antibody and the same antibody preincubated with blocking peptide.



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using BAD(Ab-112) Antibody (left) or the same antibody preincubated with blocking peptide (right).

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.