





Bad (Ab-136) Antibody

Product Code	CSB-PA810606
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q61337
Immunogen	Peptide sequence around aa.134~138 (S-R-S-A-P) derived from Mouse BAD.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total BAD protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:1000,IHC:1:50-1:200
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	Bbc2
Product Type	Polyclonal Antibody
Immunogen Species	Mus musculus (Mouse)
Target Names	Bad

Image



Western blot analysis of extracts from HT29 and HepG2 cells using BAD(Ab-136) Antibody.



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using BAD(Ab-136) Antibody(left) or the same antibody preincubated with blocking peptide(right).