




TP73 (Ab-99) Antibody

Product Code	CSB-PA938470
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15350
Immunogen	Peptide sequence around aa.97~101 (S-P-Y-A-Q) derived from Human P73.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total p73 protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	TP73; TRP73; p53-related protein; p73-alpha;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TP73

Image



Western blot analysis of extracts from HT29 and 293 cells using p73(Ab-99) Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.