





# PTK2 (Ab-861) Antibody

<b>Product Code</b>	CSB-PA939147
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q05397
<b>Immunogen</b>	Peptide sequence around aa. 859~863 (H-I-Y-Q-P) derived from Human FAK.
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Specificity</b>	The antibody detects endogenous level of total FAK protein.
<b>Tested Applications</b>	ELISA,WB,IHC;WB:1:500-1:1000,IHC:1:50-1:200
<b>Form</b>	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
<b>Purification Method</b>	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
<b>Clonality</b>	Polyclonal
<b>Alias</b>	FADK 1; FAK1; PTK2
<b>Product Type</b>	Polyclonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	PTK2

## Image



Western blot analysis of extracts from HT29, HeLa and 293 cells using FAK(Ab-861) Antibody.



Immunohistochemical analysis of paraffin-embedded human lung carcinoma tissue using FAK(Ab-861) Antibody(left) or the same antibody preincubated with blocking peptide(right).