





ITGB3 (Ab-773) Antibody

Product Code	CSB-PA243664
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05106
Immunogen	Peptide sequence around aa. 771~775 (P-L-Y-K-E) derived from Human Integrin b3.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous level of total Integrin b3 protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:1000,IHC:1:50-1:200
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	Platelet membrane glycoprotein IIIa
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ITGB3

Image



Western blot analysis of extracts from HL60 and HT29 cells using Integrin b3(Ab-773) Antibody.



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using Integrin b3(Ab-773) Antibody(left) or the same antibody preincubated with blocking peptide(right).