




TP53 (Ab-9) Antibody

Product Code	CSB-PA599282
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04637
Immunogen	Peptide sequence around aa.7~11 (D-P-S-V-E) derived from Human p53.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous level of total p53 protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	Tumor suppressor p53; Phosphoprotein p53; Antigen NY-CO-13; TP53;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TP53

Image



Western blot analysis of extracts from 293, HT29 and MDA cells using p53(Ab-9) Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.