




EZR (Ab-353) Antibody

Product Code	CSB-PA940826
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15311
Immunogen	Peptide sequence around aa.351~355 (Q-D-Y-E-E) derived from Human EZRIN.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous level of total Ezrin protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	Cytovillin; EZRI; VIL2; Villin 2; p81
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	EZR

Image



Western blot analysis of extracts from MDA and SKOV3 cells using Ezrin(Ab-353) Antibody.