




CTNNB1 (Ab-41/45) Antibody

Product Code	CSB-PA084843
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35222
Immunogen	Peptide sequence around aa.39~43/43~47 (A-T-T-T-A-P-S-L-S) derived from Human b-Catenin.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total b-Catenin protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	CTNNB1; CATNB; CTNB1; CTNNB;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CTNNB1

Image



Western blot analysis of extracts from HeLa and HT29 cells using b-Catenin(Ab-41/45) Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.