




AR (Ab-650) Antibody

Product Code	CSB-PA227367
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10275
Immunogen	Peptide sequence around aa.648~652 (T-T-S-P-T) derived from Human Androgen Receptor.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous level of total Androgen Receptor protein.
Tested Applications	ELISA, WB; WB: 1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	ANDR; DHTR; AR
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	AR

Image



Western blot analysis of extracts from Mouse kidney and brain tissue using Androgen Receptor (Ab-650) Antibody.