




ICAM1 (Ab-512) Antibody

Product Code	CSB-PA918533
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05362
Immunogen	Peptide sequence around aa. 510~514 (K-K-Y-R-L) derived from Human ICAM-1.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous level of total ICAM-1 protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	ICA1; ICAM1;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ICAM1

Image



Western blot analysis of extracts from MDA cells using ICAM-1(Ab-512) Antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.