




RB1 (Ab-795) Antibody

Product Code	CSB-PA990932
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06400
Immunogen	Peptide sequence around aa.793~797 (P-S-S-P-L) derived from Human Rb.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total Rb protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	P105-RB; PP105; PP110; RB-1; RB1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RB1

Image



Western blot analysis of extracts from MOLT cells using Rb(Ab-795) Antibody and the same antibody preincubated with blocking peptide.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.