




PRKD1 (Ab-910) Antibody

Product Code	CSB-PA989005
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15139
Immunogen	Peptide sequence around aa. 908~912 (R-V-S-I-L) derived from Human PKD/PKCm.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total PKD/PKCm protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	nPKC-D1; Protein kinase D; Protein kinase C mu type; nPKC-mu;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PRKD1

Image



Western blot analysis of extracts from HeLa, 293 and 3T3 cells using PKD/PKCμ (Ab-910) Antibody.