





# MAP3K5 (Ab-966) Antibody

<b>Product Code</b>	CSB-PA085364
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q99683
<b>Immunogen</b>	Peptide sequence around aa.964~968 (S-I-S-L-P) derived from Human ASK1.
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse
<b>Specificity</b>	The antibody detects endogenous level of total ASK1 protein.
<b>Tested Applications</b>	ELISA,WB,IHC,IF;WB:1:500-1:1000,IHC:1:50-1:200,IF:1:100-1:200
<b>Form</b>	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
<b>Purification Method</b>	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
<b>Clonality</b>	Polyclonal
<b>Alias</b>	ASK-1; M3K5; MAP3K5; MAPK/ERK kinase kinase 5; MAPKKK5
<b>Product Type</b>	Polyclonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	MAP3K5


## Image



Western blot analysis of extracts from 293 and HeLa cells using ASK1 (Ab-966) Antibody.



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using ASK1 (Ab-966) Antibody (left) or the same antibody preincubated with blocking peptide (right).



Immunofluorescence staining of methanol-fixed HeLa cells using ASK1 (Ab-966) Antibody.

---

**Usage**

For Research Use Only. Not for use in diagnostic or therapeutic procedures.