





Histone H3.1 (Ab-10) Antibody

Product Code	CSB-PA585914
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P68431
Immunogen	Peptide sequence around aa.8~12 (R-K-S-T-G) derived from Human Histone H3.1.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total Histone H3.1 protein.
Tested Applications	ELISA,WB,IF;WB:1:500-1:1000,IF:1:100-1:200
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	H3/a; H3/c; H3/d; H3/f; H3/h
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	HIST1H3A

Image



Western blot analysis of extracts from 293 cells using Histone H3.1(Ab-10) Antibody.



Immunofluorescence staining of methanol-fixed HeLa cells using Histone H3.1(Ab-10) Antibody.



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.