





HDAC4 (Ab-632) Antibody

Product Code	CSB-PA078369
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P56524
Immunogen	Peptide sequence around aa. 630~634 (A-Q-S-S-P) derived from Human HDAC4.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total HDAC4 protein.
Tested Applications	ELISA,WB,IF;WB:1:500-1:1000,IF:1:100-1:200
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	HD4; Histone deacetylase 4;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	HDAC4

Image



Western blot analysis of extracts from HepG2 cells using HDAC4(Ab-632) Antibody.



Immunofluorescence staining of methanol-fixed HeLa cells using HDAC4(Ab-632) Antibody.



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.