





CDKN1B (Ab-10) Antibody

Product Code	CSB-PA018782
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P46527
Immunogen	Peptide sequence around aa.8~12 (N-G-S-P-S) derived from Human p27Kip1.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total p27Kip1 protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:1000,IHC:1:50-1:200
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	CDKN1B; CDN1B; Cyclin-dependent kinase inhibitor 1B; Cyclin-dependent kinase inhibitor p27; KIP1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CDKN1B

Image



Western blot analysis of extracts from A2780 cells using p27Kip1(Ab-10)



Immunohistochemical analysis of paraffin-embedded human lung carcinoma tissue using p27Kip1(Ab-10)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.