




GRIN2B (Ab-1474) Antibody

Product Code	CSB-PA113700
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13224
Immunogen	Peptide sequence around aa.1472~1476 (H-V-Y-E-K) derived from Human NMDAR2B.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total NMDAR2B protein.
Tested Applications	ELISA,IF;IF:1:100-1:200
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	GRIN2B; NMDE2; NME2; NR2B; NR3
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GRIN2B

Image



Immunofluorescence staining of methanol-fixed HeLa cells using NMDAR2B(Ab-1474) antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.