





NOS1 (Ab-852) Antibody

Product Code	CSB-PA149037
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29475
Immunogen	Peptide sequence around aa.850~854 (F-N-S-V-S) derived from Human nNOS.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total nNOS protein.
Tested Applications	ELISA,WB,IF;WB:1:500-1:1000,IF:1:100-1:200
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	BNOS; Constitutive NOS; N-NOS; NC-NOS; NOS1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NOS1

Image



Western blot analysis of extracts from Rat and Mouse brain tissue using nNOS(Ab-852) Antibody.



Immunofluorescence staining of methanol-fixed HeLa cells using nNOS(Ab-852) Antibody.

**Usage**

For Research Use Only. Not for use in diagnostic or therapeutic procedures.