





FOXO4 (Ab-197) Antibody

Product Code	CSB-PA997450
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P98177
Immunogen	Peptide sequence around aa. 195~199 (A-A-S-M-D) derived from Human AFX.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total AFX protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:1000,IHC:1:50-1:200
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	AFX; FOXO4; AFX1; Afxh;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FOXO4

Image



Western blot analysis of extracts from 293 and THP1 cells using AFX(Ab-197) Antibody.



Immunohistochemical analysis of paraffin-embedded human lung carcinoma tissue using AFX(Ab-197) Antibody(left) or the same antibody preincubated with blocking peptide(right).