




VASP (Ab-239) Antibody

Product Code	CSB-PA206698
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50552
Immunogen	Peptide sequence around aa. 236~240 (K-V-S-K-Q) derived from Human VASP.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous level of total VASP protein.
Tested Applications	ELISA, WB; WB: 1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	VASP

Image



Western blot analysis of extracts from 293 and C6 cells using VASP (Ab-239) Antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.