




PRKDC (Ab-2609) Antibody

Product Code	CSB-PA178115
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P78527
Immunogen	Peptide sequence around aa.2607~2611 (V-E-T-Q-A) derived from Human DNA-PK.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous level of total DNA PKcs protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	DNPK1; PRKD; PRKDC; XRCC7; P460
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PRKDC

Image



Western blot analysis of extracts from HT29 and HepG2 cells using DNA PKcs(Ab-2609) antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.