





# FLT3 (Ab-591) Antibody

<b>Product Code</b>	CSB-PA546306
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P36888
<b>Immunogen</b>	Peptide sequence around aa.589~593 (Y-F-Y-V-D) derived from Human FLT3.
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Specificity</b>	The antibody detects endogenous level of total FLT3 protein.
<b>Tested Applications</b>	ELISA,WB,IF;WB:1:500-1:1000,IF:1:100-1:200
<b>Form</b>	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
<b>Purification Method</b>	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
<b>Clonality</b>	Polyclonal
<b>Alias</b>	CD135; FLK-2; FLT-3; FLT3; STK-1
<b>Product Type</b>	Polyclonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	FLT3

## Image



Western blot analysis of extracts from 293 and MCF cells using FLT3(Ab-591) Antibody.



Immunofluorescence staining of methanol-fixed Hela cells using FLT3(Ab-591) Antibody.

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.