





FLT3 (Ab-591) Antibody

Product Code	CSB-PA546306
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P36888
Immunogen	Peptide sequence around aa.589~593 (Y-F-Y-V-D) derived from Human FLT3.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous level of total FLT3 protein.
Tested Applications	ELISA, WB, IF; WB: 1:500-1:1000, IF: 1:100-1:200
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	CD135; FLK-2; FLT-3; FLT3; STK-1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FLT3

Image



Western blot analysis of extracts from 293 and MCF cells using FLT3(Ab-591) Antibody.



Immunofluorescence staining of methanol-fixed HeLa cells using FLT3(Ab-591) Antibody.