




EFNB2 (Ab-316) Antibody

Product Code	CSB-PA587215
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P52799
Immunogen	Peptide sequence around aa.314~318 (P-V-Y-I-V) derived from Human Ephrin B (ephrin-B2).
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous level of total Ephrin-B2 protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	EFNB2; HTKL; LERK5
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	EFNB2

Image



Western blot analysis of extracts from HUVEC cells using Ephrin-B2(Ab-316) Antibody and the same antibody preincubated with blocking peptide.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.