





EFNB2 (Ab-330) Antibody

Product Code	CSB-PA077593
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P52799
Immunogen	Peptide sequence around aa.328~332 (N-I-Y-Y-K) derived from Human Ephrin B (ephrin-B2).
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous level of total Ephrin-B2 protein.
Tested Applications	ELISA,WB,IF;WB:1:500-1:1000,IF:1:100-1:200
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	EFNB2; HTKL; LERK5
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	EFNB2

Image



Western blot analysis of extracts from HT29 cells using Ephrin-B2(Ab-330) Antibody.



Immunofluorescence staining of methanol-fixed HeLa cells using Ephrin-B2(Ab-330) Antibody.