





# RAC1 (Ab-71) Antibody

<b>Product Code</b>	CSB-PA289719
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P63000
<b>Immunogen</b>	Peptide sequence around aa.69~73 (P-L-S-Y-P) derived from Human Rac1.
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Specificity</b>	The antibody detects endogenous level of total Rac1 protein.
<b>Tested Applications</b>	ELISA,WB,IF;WB:1:500-1:1000,IF:1:100-1:200
<b>Form</b>	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
<b>Purification Method</b>	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Ras-like protein TC25; Ras-related C3 botulinum toxin substrate 1; p21-Rac1
<b>Product Type</b>	Polyclonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	RAC1

## Image



Western blot analysis of extracts from 293 and HeLa cells using Rac1(Ab-71) Antibody .



Immunofluorescence staining of methanol-fixed HeLa cells using Rac1(Ab-71) Antibody.