





RAC1 (Ab-71) Antibody

Product Code	CSB-PA289719
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P63000
Immunogen	Peptide sequence around aa.69~73 (P-L-S-Y-P) derived from Human Rac1.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total Rac1 protein.
Tested Applications	ELISA,WB,IF;WB:1:500-1:1000,IF:1:100-1:200
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	Ras-like protein TC25; Ras-related C3 botulinum toxin substrate 1; p21-Rac1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RAC1

Image



Western blot analysis of extracts from 293 and HeLa cells using Rac1(Ab-71) Antibody .



Immunofluorescence staining of methanol-fixed HeLa cells using Rac1(Ab-71) Antibody.



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.