





APP (Ab-668) Antibody

Product Code	CSB-PA907645
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05067
Immunogen	Peptide sequence around aa.666~670 (A-V-T-P-E) derived from Human APP.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total APP protein.
Tested Applications	ELISA,WB,IF;WB:1:500-1:1000,IF:1:100-1:200
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	AAA; AD1; PN2; ABPP; APPI
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	APP

Image



Western blot analysis of extracts from mouse brain tissue using APP(Ab-668) Antibody and the same antibody preincubated with blocking peptide.



Immunofluorescence staining of methanol-fixed MCF cells using APP(Ab-668) Antibody.



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.