




GABRB1 (Ab-434) Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA620050 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P18505 |
| Immunogen | Peptide sequence around aa.432~436 (R-A-S-Q-L) derived from Human GABRB1. |
| Raised In | Rabbit |
| Species Reactivity | Human,Mouse,Rat |
| Specificity | The antibody detects endogenous level of total GABRB1 protein. |
| Tested Applications | ELISA,WB;WB:1:500-1:1000 |
| Form | Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Purification Method | Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide. |
| Clonality | Polyclonal |
| Alias | GABA(A) receptor; GABA(A)R beta-1; GBRB1; GABA-RB; |
| Product Type | Polyclonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | GABRB1 |

Image



Western blot analysis of extracts from PC12 cells using GABRB1(Ab-434) Antibody and the same antibody preincubated with blocking peptide.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.