




GABRB1 (Ab-434) Antibody

Product Code	CSB-PA620050
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P18505
Immunogen	Peptide sequence around aa.432~436 (R-A-S-Q-L) derived from Human GABRB1.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total GABRB1 protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	GABA(A) receptor; GABA(A)R beta-1; GBRB1; GABA-RB;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GABRB1

Image



Western blot analysis of extracts from PC12 cells using GABRB1(Ab-434) Antibody and the same antibody preincubated with blocking peptide.