





VASP (Ab-157) Antibody

Product Code	CSB-PA158569
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50552
Immunogen	Peptide sequence around aa.155~159 (R-V-S-N-A) derived from Human VASP.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous level of total VASP protein.
Tested Applications	ELISA,WB,IF;WB:1:500-1:1000,IF:1:100-1:200
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	VASP

Image



Western blot analysis of extracts from 293 and C6 cells using VASP(Ab-157) Antibody.



Immunofluorescence staining of methanol-fixed HeLa cells using VASP(Ab-157) Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.