




# Met Antibody

<b>Product Code</b>	CSB-PA047827
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P08581
<b>Immunogen</b>	Peptide sequence around aa.1347~1351 (E-H-Y-V-H) derived from Human Met.
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Tested Applications</b>	ELISA,WB;WB:1:500-1:1000
<b>Form</b>	Liquid
<b>Purification Method</b>	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
<b>Clonality</b>	Polyclonal
<b>Product Type</b>	Polyclonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	MET

## Image



Western blot analysis of extracts from HepG2 cells using Met(Ab-1349) Antibody and the same antibody preincubated with blocking peptide.