




Met Antibody

Product Code	CSB-PA047827
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08581
Immunogen	Peptide sequence around aa.1347~1351 (E-H-Y-V-H) derived from Human Met.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Liquid
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MET

Image



Western blot analysis of extracts from HepG2 cells using Met(Ab-1349) Antibody and the same antibody preincubated with blocking peptide.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.