





BRCA1 (Ab-1423) Antibody

Product Code	CSB-PA186422
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P38398
Immunogen	Peptide sequence around aa.1421~1425 (H-G-S-Q-P) derived from Human BRCA1.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous level of total BRCA1 protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:1000,IHC:1:50-1:200
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	RNF53
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	BRCA1

Image



Western blot analysis of extracts from HeLa and HT29 cells using BRCA1(Ab-1423) Antibody.



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using BRCA1(Ab-1423) Antibody(left) or the same antibody preincubated with blocking peptide(right).

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.