




NOS3 (Ab-495) Antibody

Product Code	CSB-PA064019
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29474
Immunogen	Peptide sequence around aa.493~497 (K-K-T-F-K) derived from Human eNOS.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total eNOS protein.
Tested Applications	ELISA,IF;IF:1:100-1:200
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	Constitutive NOS; EC-NOS; ECNOS; NOS3; NOSIII
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NOS3

Image



Immunofluorescence staining of methanol-fixed HeLa cells using eNOS(Ab-495) Antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.