




# NOS3 (Ab-495) Antibody

<b>Product Code</b>	CSB-PA064019
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P29474
<b>Immunogen</b>	Peptide sequence around aa.493~497 (K-K-T-F-K) derived from Human eNOS.
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Specificity</b>	The antibody detects endogenous level of total eNOS protein.
<b>Tested Applications</b>	ELISA,IF;IF:1:100-1:200
<b>Form</b>	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
<b>Purification Method</b>	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Constitutive NOS; EC-NOS; ECNOS; NOS3; NOSIII
<b>Product Type</b>	Polyclonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	NOS3

## Image



Immunofluorescence staining of methanol-fixed HeLa cells using eNOS(Ab-495) Antibody.

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.