




MAPK14 (Ab-182) Antibody

Product Code	CSB-PA112086
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16539
Immunogen	Peptide sequence around aa. 180~184 (T-G-Y-V-A) derived from Human P38 MAPK.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total P38MAPK protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	MAPK2; MAPKAPK-2; MAPKAPK2
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MAPK14

Image



Western blot analysis of extracts from 293 and HeLa cells using P38 MAPK(Ab-182) Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.